




Rapid Rhino®



A Full Range for Epistaxis

Easy to insert and remove due to self-lubricating hydrocolloid fabric and ultra-low profile.



RR 450 Model



Anterior

RR 450 4.5cm

RR 550 5.5cm

RR 750 7.5cm

Anterior with Airway

RR 551 5.5cm

RR 751 7.5cm



The central airway of the **RR 551 & RR 751** allows the patient some continued nasal airflow with the device properly positioned.

ANTERIOR/ANTERIOR WITH AIRWAY

Quickly conforms to nasal anatomy and provides gentle and even compression to areas of epistaxis.



RR 552 Model



Bilateral

| RR 552 5.5cm

| RR 752 7.5cm

Bilateral with Airway

| RR 555 5.5cm

| RR 755 7.5cm




The central airway of the **RR 555 & RR 755** allows the patient some continued nasal airflow with the device properly positioned.


BILATERAL/BILATERAL WITH AIRWAY

Contains two balloons to control intranasal hemorrhage originating from the sphenopalatine artery for posterior epistaxis.

Inflating pilot cuffs provide tactile feel to easily gauge balloon pressure during inflation

CMC-infused fabric exterior wets into a self-lubricating gel for optimal patient comfort from insertion to removal





AIR inflation provides a delicate and lightweight pressure mechanism

Green-lined swallow guard provides indication for posterior balloon

Two independent balloons individually conform to the anatomy and create direct pressure to bleeds

Posterior

| RR 900 9.0cm

POSTERIOR

A Full Range for Epistaxis.

Be prepared for your epistaxis procedures with Rapid Rhino[®], the product with the widest range of sizes for the widest range of epistaxis events.

Stock Rapid Rhino and you will always be ready.





ENT REP BC
FPO

Check out our full range of Rapid Rhino epistaxis products at:
www.rapidrhino.com or call **1.800.797.6520** for more information.



ArthroCare ENT

680 Vaqueros Avenue
Sunnyvale, CA 94085-3523
main phone 408-736-0224
order entry **800-797-6520**
order entry fax 888-994-2782
www.arthrocareENT.com

ArthroCare Europe AB

Skeppsbron 2
111 30 Stockholm
Sweden
phone +46 8 546 172 00
fax +46 8 546 172 39
info@arthrocare.se

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

©2010 ArthroCare Corporation. All rights reserved. All trademarks, service marks, and logos listed herein are registered and/or unregistered trademarks of ArthroCare Corporation, its affiliates, subsidiaries, or a third party who has licensed its trademarks to ArthroCare. ArthroCare ENT, ArthroCare Sports Medicine, and ArthroCare Spine are affiliates of ArthroCare Medical Corporation.

▶ **Stock Rapid Rhino**
and you will always
be ready.

